



Belegung	FS-CW-624 <i>ID 0,32mm, df 1,80µm, Länge 60m</i> 6% Cyanopropylsiloxan Max. Temperatur: 240/260°C
Dimensionen	60 m x 0,32 mm ID, df 1,8 µm
Artikel Nr.	23214360
Probe	2,3-Dimethyl-1-buten, 2,3-Dimethylbutan, Cyclohexan
Ofen	50 °C, 10 min, 5 °C/min, 120 °C
Trärgas	1 bar N2
Detektor	FID, 320 °C

