



Belegung

FS-SE-54

Dimensionen

2 x 50 m x 0,32 mm ID (gekoppelt)

Artikel Nr.

[232057100](#)

Probe

Trennung Alkene (C2-C4)

Substanzen: Ethen / Propen / trans-Buten / iso-Buten / cis-Buten

Ofen

60°C isotherm

Trärgas

0,8 bar N₂, Split geschlossen

Detektor

FID 200°C, 1 x 16



1.) Ethen

4.) Isobuten

2.) Propen

5.) cis-Buten

3.) trans-Buten