



Belegung	FS-OV-1701 <i>ID 0,25mm, df 0,25µm, Länge 50m</i> <i>86% Methyl- 7% cyanopropyl- 7%</i> <i>phenylsilikon</i> <i>Max. Temperatur: 260°C</i>
Dimensionen	50 m x 0,25 mm ID, df 0,25 µm
Artikel Nr.	<a href="#">225060</a>
Probe	BTX-Aromaten in Dichlormethan
Ofen	80°C isotherm
Trägergas	1,0 bar N2
Detektor	FID 250°C

