



Belegung

FS-SE-54

ID 0,25mm, df 0,50µm, Länge 50m
94% Methyl- 5% phenyl- 1% vinylsilikon
Max. Temperatur: 300/320°C

Dimensionen

50 m x 0,25 mm ID, df 0,5 µm

Artikel Nr.

[225054](#)

Probe

PCB in A40/A60


Substanzen: PCB 31 / PCB 28 / PCB 52 / PCB 101 / PCB 153 / PCB 138 / PCB 180

Ofen

Trärgas

Detektor

ECD



1.) PCB's