




Belegung	FS-SE-54 <i>ID 0,32mm, df 1µm, Länge 25m</i> <i>94% Methyl- 5% phenyl- 1% vinylsilikon</i> <i>Max. Temperatur: 300/320°C</i>
Dimensionen	25 m x 0,32 mm ID, df 1,21 µm
Artikel Nr.	<a href="#">232056</a>
Probe	Dioxan in Dichlormethan
Ofen	80°C
Trägergas	0,7 bar N <sub>2</sub>
Detektor	FID 250°C



1.) Dichlormethan

2.) Dioxan